

BEYOND MOTION



# CAPSULE

CAMERA ENCLOSURE

NICON

# THE WEATHERPROOF SOLUTION FOR OUTDOOR MOTION CAPTURE

Developed for use in a large-scale permanent desert installation, Vicon Capsule is also ideal for university researchers, biomechanics tracking in large sports arenas, automotive companies looking to capture motion on outdoor tracks and for aerospace applications.

Independently verified to IP65 rating and weatherproof, Capsule protects Vantage cameras from even the most adverse weather conditions.

Dust-tight and water resistant, it protects your Vantage camera in temperatures from -5°C to 35°C and can be left permanently installed outdoors, exposed to the elements, allowing for continual capture, without delay.

Capsule supports all Vantage camera models, as well as ViperX too.


Capsule also reduces the unit's weight by over 50 percent to 1.7kg (3.3kg including the camera) and introduces new features for safer carrying and installation, including tether fittings and a handle.

You can capture the same volume with precisely the same setup for regular, repeat testing over days, weeks, months and even years.

It enables you to react to changing capture needs within complex projects, without having to stop to set up the system. With Capsule, the Vicon system is ready when you are.



CAPSULE  
CAMERA ENCLOSURE

CAPSULE		Technical Specifications	
 <b>VICON</b> Contact us today to book a demo or find out more. <a href="mailto:info@vicon.com">info@vicon.com</a> <a href="http://vicon.com">vicon.com</a> © 2021 Vicon Motion Systems Ltd. All rights reserved. UK Registered no. 1801446	RATING	IP65	
	WEIGHT	1.7 kg (3.3 kg inc camera)	
	DIMENSIONS	140 x 180 x 200 mm	
	CAMERAS SUPPORTED	ViperX, VANTAGE V5, V8, V16 with variable lenses	